



# SILICONE PENROSE DRAIN TUBE




# SILICONE PENROSE DRAIN TUBE

Fortune® Silicone Penrose Drainage Tube are used for post-operative drainage. Soft and smooth tube minimizes trauma during insertion and withdrawal, the tube draws the fluids to the exterior using the capillary phenomenon. Transparent tube with X-ray opaque line allows for easy visual inspection, fluid observation, and also confirmation of intubated tube using X-ray; white tube has X-ray opaque effectiveness throughout the entire length. Secure-Clip provides the tube an easy fixation on the wound site.

Three types of Penrose are available for selection:

- Penrose Drain Tube
- Penrose Drain Tube w/ Secure-Clip
- Penrose Straight Drain Tube

- A** Silicone Penrose Drain Tube
- B** Silicone Penrose Drain Tube w/ Secure-Clip
- C** Secure-Clip
- D** Silicone Penrose Straight Drain Tube




## SPECIFICATIONS

### Penrose Drain Tube w/ Secure-Clip:

REF. No.	O.D.	Description
2010-5405	5 mm	
2010-5406	6 mm	
2010-5408	8 mm	- WS-type - 300 mm
2010-5410	10 mm	- White opaque effectiveness
2010-5412	12 mm	- Secure-clip
2010-5415	15 mm	- Graduation
2010-5419	19 mm	
2010-5425	25 mm	

REF. No.	O.D.	Description
2010-6405	5 mm	
2010-6406	6 mm	
2010-6408	8 mm	- WS-type - 500 mm
2010-6410	10 mm	- White opaque effectiveness
2010-6412	12 mm	- Secure-clip
2010-6415	15 mm	- Graduation
2010-6419	19 mm	
2010-6425	25 mm	



2010-54XX/64XX


# SILICONE PENROSE DRAIN TUBE

### Penrose Drain Tube:

REF. No.	O.D.	Description
2010-0405	5 mm	
2010-0406	6 mm	
2010-0408	8 mm	- S-type
2010-0410	10 mm	- 300 mm
2010-0412	12 mm	- Transparent
2010-0412	12 mm	- X-ray opaque line
2010-0415	15 mm	
2010-0419	19 mm	
2010-0425	25 mm	
2010-1405	5 mm	
2010-1406	6 mm	
2010-1408	8 mm	- S-type
2010-1410	10 mm	- 500 mm
2010-1410	10 mm	- Transparent
2010-1412	12 mm	- X-ray opaque line
2010-1415	15 mm	
2010-1419	19 mm	
2010-1425	25 mm	

### Penrose Straight Drain Tube:

REF. No.	Size	Width	Length	Description
2011-0020	S	20 mm		
2011-0025	M	25 mm	400 mm	X-ray opaque effectiveness
2011-0035	L	35 mm		
2011-0045	XL	45 mm		



### Secure-Clip:

REF. No.	Size	Description
2010-0000	Universal	Metal



2010-04XX/14XX



2010-0000



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,  
 Danshuei Dist, New Taipei City. 251 Taiwan  
 Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266  
 www.fortunemed.com Mail:info@fortunemed.com