



SILICONE T DRAINAGE TUBE

Choices for alternative drainage



GENERAL
SURGICAL
SERIES




SILICONE T DRAINAGE TUBE

Fortune® Silicone T Drainage Tube Feature:

- Medical grade silicone for superior biocompatibility.
- Proper rigidity design to prevent from kink or clog over time.
- Larger diameter of connecting tube allows smooth withdrawn.
- X-ray opaque line.
- Sterilized double packaging.

A Silicone T Drainage Tube

B Silicone T Drainage Tube Set




SPECIFICATIONS

T Drainage Tube:


REF. No.	Size	Description
1910-0005	5 FR	
1910-0008	8 FR	
1910-0010	10 FR	
1910-0012	12 FR	- A: 180 mm
1910-0014	14 FR	- B: 450 mm
1910-0016	16 FR	- X-ray opaque line
1910-0018	18 FR	- w/large bore male connector
1910-0020	20 FR	
1910-0022	22 FR	
1910-0024	24 FR	

T Drainage Tube With Eyes:

REF. No.	Size	Description
1910-0210	10 FR	
1910-0212	12 FR	
1910-0214	14 FR	- Perforation
1910-0216	16 FR	- A: 180 mm
1910-0218	18 FR	- B: 450 mm
1910-0220	20 FR	- X-ray opaque line
1910-0222	22 FR	- w/large bore male connector
1910-0224	24 FR	



1910-00XX



1910-02XX


SILICONE T DRAINAGE TUBE

T Drainage Set:


REF. No.	Size	Description
1910-4005	5 FR	
1910-4008	8 FR	- T Drainage Tube
		- Connector
1910-4010	10 FR	- 400 cc Collection Bag
1910-8012	12 FR	
1910-8014	14 FR	
1910-8016	16 FR	
1910-8018	18 FR	
		- T Drainage Tube
1910-8020	20 FR	- Connector
		- 800 cc Collection Bag
1910-8022	22 FR	
1910-8024	24 FR	

T Drainage Set With Eyes:

REF. No.	Size	Description
1910-8212	12 FR	
1910-8214	14 FR	
1910-8216	16 FR	- T Drainage Tube With Eyes
		- Connector
		- 800 cc Collection Bag
1910-8218	18 FR	
1910-8220	20 FR	
1910-8222	22 FR	
1910-8224	24 FR	



1910-80XX



1910-82XX



fortune
medical

6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com

