



# SILICONE HEMATURIA CATHETER

Your Right Choice for Quality Catheters



UROLOGICAL  
SERIES





# SILICONE HEMATURIA CATHETER

## Fortune® Silicone 3-Way Hematuria Catheter Features:

- Medical grade silicone for superior biocompatibility.
- Transparent tube allows easy visual inspection and fluid observation.
- X-ray opaque line allows intubating location confirmation.
- A soft and uniformly inflated larger volume balloon makes the tube sit well against the bladder.
- 3-Way design with an attached plug on the irrigation funnel.
- Tapered Central funnel could be connected with irrigation bags.
- Five types are available for various applications:
  - A. type: Straight round tip
  - B. type: Straight whistle tip
  - C. type: Coude round tip
  - D. type: Dufour
  - E. type: Couvelaire
- Sterilized double packaging.


- A** Silicone Hematuria Catheter, Straight Round Tip
- B** Silicone Hematuria Catheter, Straight Whistle Tip
- C** Silicone Hematuria Catheter, Coude Round Tip
- D** Silicone Hematuria Catheter, Dufour
- E** Silicone Hematuria Catheter, Couvelaire



## SPECIFICATIONS

### Hematuria Catheter, Straight Round Tip:


REF. No.	Size	Balloon	Length	Description
1861-0318	18 FR	50-80 cc/ml		- Straight round tip
1861-0320	20 FR		420 mm	- X-ray opaque line
1861-0322	22 FR	80-100 cc/ml		- Attached plug
1861-0324	24 FR			



# SILICONE HEMATURIA CATHETER

## Hematuria Catheter, Straight Whistle Tip:


REF. No.	Size	Balloon	Length	Description
1862-0318	18 FR	50-80 cc/ml		- Straight whistle tip
1862-0320	20 FR		420 mm	- X-ray opaque line
1862-0322	22 FR	80-100 cc/ml		- Attached plug
1862-0324	24 FR			



1862-03XX

## Hematuria Catheter, Coude Round Tip:


REF. No.	Size	Balloon	Length	Description
1863-0318	18 FR	50-80 cc/ml		- Coude round tip
1863-0320	20 FR		420 mm	- X-ray opaque line
1863-0322	22 FR	80-100 cc/ml		- Attached plug
1863-0324	24 FR			



1863-03XX

## Hematuria Catheter, Dufour:


REF. No.	Size	Balloon	Length	Description
1864-0318	18 FR	50-80 cc/ml		- Dufour
1864-0320	20 FR		420 mm	- X-ray opaque line
1864-0322	22 FR	80-100 cc/ml		- Attached plug
1864-0324	24 FR			



1864-03XX

## Hematuria Catheter, Couvelaire:

REF. No.	Size	Balloon	Length	Description
1865-0318	18 FR	50-80 cc/ml		- Couvelaire
1865-0320	20 FR		420 mm	- X-ray opaque line
1865-0322	22 FR	80-100 cc/ml		- Attached plug
1865-0324	24 FR			



1865-03XX



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,  
 Danshuei Dist, New Taipei City, 251 Taiwan  
 Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266  
 www.fortunemed.com Mail:info@fortunemed.com