



# SILICONE SUPRA-PUBIC CATHETER

Quality, Safety & Reliability Guaranteed




UROLOGICAL  
SERIES

# SILICONE SUPRA-PUBIC CATHETER

## Fortune® Silicone Supra-Pubic Catheter Features:

- Medical grade silicone for superior biocompatibility.
- Transparent with X-ray opaque line.
- Integrated balloon for smooth insertion through split cannula or surgical procedure.
- Smooth round shaft minimize trauma during insertion and withdrawal.
- Pigtail pattern available for enhanced drainage.
- Sterilized double packaging.

- A Silicone Supra-pubic Catheter**
- B Silicone Supra-Pubic Catheter Puncture Set**
- C Silicone Supra-Pubic Catheter Exchange Set**



## SPECIFICATIONS

### Supra-pubic Catheter:


| Ordinary Balloon REF. No. | Integrated Balloon REF. No. | Size  | Balloon    | Length | Description     |
|---------------------------|-----------------------------|-------|------------|--------|-----------------|
| 1890-0212                 | 4890-0212                   | 12 FR |            |        |                 |
| 1890-0214                 | 4890-0214                   | 14 FR |            |        |                 |
| 1890-0216                 | 4890-0216                   | 16 FR |            |        |                 |
| 1890-0218                 | 4890-0218                   | 18 FR | 5-10 cc/ml | 420 mm | - Depth marking |
| 1890-0220                 | 4890-0220                   | 20 FR |            |        |                 |
| 1890-0222                 | 4890-0222                   | 22 FR |            |        |                 |
| 1890-0224                 | 4890-0224                   | 24 FR |            |        |                 |



1890-02XX

### Supra-Pubic Catheter Puncture Set:

| REF. No.  | Size  | Balloon    | Length | Description  |
|-----------|-------|------------|--------|--|
| 4890-9212 | 12 FR | 5-10 cc/ml | 420 mm | - Integrated balloon<br>- Split cannula<br>- Depth marking |
| 4890-9214 | 14 FR |            |        | - Catheter clamp   |




4890-92XX

# SILICONE SUPRA-PUBIC CATHETER

**Supra-Pubic Catheter Exchange Set:**

| Ordinary Balloon REF. No. | Integrated Balloon REF. No. | Size  | Balloon    | Length | Description                    |
|---------------------------|-----------------------------|-------|------------|--------|--------------------------------|
| 1890-8212                 | 4890-8212                   | 12 FR |            |        |                                |
| 1890-8214                 | 4890-8214                   | 14 FR |            |        |                                |
| 1890-8216                 | 4890-8216                   | 16 FR |            |        | - Depth marking                |
| 1890-8218                 | 4890-8218                   | 18 FR | 5-10 cc/ml | 420 mm | - Catheter clamp               |
| 1890-8220                 | 4890-8220                   | 20 FR |            |        | - PTFE coated guidewire 0.035" |
| 1890-8222                 | 4890-8222                   | 22 FR |            |        |                                |
| 1890-8224                 | 4890-8224                   | 24 FR |            |        |                                |



1890-82XX



**fortune  
medical**

6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,  
 Danshuei Dist, New Taipei City. 251 Taiwan  
 Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266  
 www.fortunemed.com Mail:info@fortunemed.com

