



FORTUNE FIXPAD ANCHORING DEVICE

Comfortable support, giving patients
the freedom to move around freely
and safely.



UROLOGICAL
SERIES




FORTUNE FIXPAD® ANCHORING DEVICE

Fortune® FixPad® Anchoring Device:

Fortune® FixPad® is designed to provide for fixing the catheter tubes and make the catheter tubes are placed orderly. The FixPad® provides comfortable support, giving patients the freedom to move around freely and safely.


Features:

- Applicable to the round catheter products, the outer diameter is between 3mm to 12mm (that is, 9FR to 36FR) or an equal circumference catheter. For example, a Foley catheter, various catheter products, etc.
- The adhesive pad is made of nonwoven fabric and coated with medical pressure sensitive adhesive.
- Easy to apply, remove and reposition. Leaves no sticky residue.
- No skin preparation required before application, and no need for solutions to dissolve the adhesive when removed.
- Ideal for repeated application over the same skin area.
- Reduced risk of allergic reactions.
- Showerproof and breathable.



SPECIFICATIONS

REF. No.	Size	Outer Diameter
1800-0005	106 x 65 mm	3mm to 12mm (9FR to 36FR)



1800-0005

SILICONE FOLEY BALLOON CATHETER WITH FIXPAD®


Fortune® Silicone 2-Way Foley Integrated Balloon Catheter with FIXPAD®

INTENDED USE:

The 2-Way Foley integrated Balloon Catheter is intended for use as a urinary catheter to pass fluids to and from the urinary bladder. Fortune® FixPad® is intended for fixing the catheter tubes, provides to reduce the chance of the catheter falling off due to the patient's activity.

Features:

- Transparent tube with X-ray opaque line, transparent for easy visual inspection and fluid observation, x-ray opaque for easy confirmation of intubated tubes using X-ray.
- The soft and uniformly inflated balloon helps the tube sit comfortably against the bladder.
- Smooth round shaft reduces the chance of trauma during insertion and withdrawal.
- Firm yet flexible tip designed to aid insertion.
- Recommended usage up to 29 days for 4822 series.



SPECIFICATIONS

Silicone 2-Way Foley Integrated Balloon Catheter Sets

REF. No.	Size	Balloon	Length	Description
4822-3512	12 FR			
4822-3514	14 FR			
4822-3516	16 FR			
4822-3518	18 FR	5-10 cc/ml	420 mm	- Integrated balloon
4822-3520	20 FR			- X-ray opaque line
4822-3522	22 FR			- FixPad® anchoring device
4822-3524	24 FR			
4822-3526	26 FR			



4822-35XX