



# CATHETER SPIGOT

Close The Foley Catheter Securely and Hygienically



UROLOGICAL  
SERIES



# CATHETER SPIGOT

### Catheter Spigot:

Catheter Spigot is used to stop flow for catheters hygienically. It is non-invasive sealing the catheter for a short time to allow urine to collect in the bladder.

### Intended Use:

For closing catheters and tubes that control urination, thereby training the activity of the sphincter (self-voiding).

### Features:

- A plastic stopper for catheters provided sterile in sealed, individual packages.
- Universal size is molded to fit all sizes of catheters or tubes.
- The structure of spigot is strong and secure.
- The spigot is a single use item and is supplied sterile.



## SPECIFICATIONS

| REF. No.  | Size      | Description                                     |
|-----------|-----------|---|
| 1820-0000 | Universal | - ABS<br>- Sterilized packaging<br>- 50 pcs/box |



1820-0000



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,  
 Danshuei Dist, New Taipei City. 251 Taiwan  
 Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266  
 www.fortunemed.com Mail:info@fortunemed.com