



SILICONE ENDOTRACHEAL TUBE

Silicone for Great Comfort and Flexibility




ANESTHETIC
SERIES

SILICONE ENDOTRACHEAL TUBE

Fortune® Silicone Endotracheal Tube is made of medical grade silicone. It provides superior biocompatibility and flexibility. Great transparency allows easy visual inspection. Graduation and X-ray opaque line help depth and location confirmation. A pilot balloon is equipped to monitor balloon status. Low pressure cuff models are available to reduce the pressure on trachea wall and to provide greater comfort to patient.

- A** Silicone Endotracheal Tube, Cuffless
- B** Silicone Endotracheal Tube, OC Type
- C** Silicone Endotracheal Tube, LC Type
- D** Silicone Endotracheal Tube, OW Type
- E** Silicone Endotracheal Tube, LW Type



SPECIFICATIONS

Endotracheal Tube, Cuffless:

PED-Magill REF. No.	Size/I.D.	O.D.	Length	Description
1510-0025	2.5 mm	4.0 mm	140 mm	
1510-0030	3.0 mm	4.7 mm	160 mm	- Magill Series: Magill tip
1510-0035	3.5 mm	5.3 mm	180 mm	- Cuffless
1510-0040	4.0 mm	6.0 mm	200 mm	- Wire reinforced
1510-0045	4.5 mm	6.7 mm	220 mm	- Graduation
1510-0050	5.0 mm	7.3 mm	240 mm	- X-ray opaque line



1510-00XX







SILICONE ENDOTRACHEAL TUBE

Endotracheal Tube, OC Type (Ordinary cuff, Curved):

OC-Magill REF. No.	OC-Murphy REF. No.	Size/I.D.	O.D.	Length	Description
1520-0045	1525-0045	4.5 mm	6.7 mm	220 mm	
1520-0050	1525-0050	5.0 mm	7.3 mm	240 mm	
1520-0055	1525-0055	5.5 mm	8.0 mm	270 mm	
1520-0060	1525-0060	6.0 mm	8.7 mm	280 mm	*OC-Magill series: Magill tip
1520-0065	1525-0065	6.5 mm	9.3 mm	290 mm	*OC-Murphy series: Murphy eye
1520-0070	1525-0070	7.0 mm	10.0 mm	300 mm	
1520-0075	1525-0075	7.5 mm	10.7 mm	310 mm	Ordinary cuff Curved Graduation Pilot balloon X-ray opaque line
1520-0080	1525-0080	8.0 mm	11.3 mm	320 mm	
1520-0085	1525-0085	8.5 mm	12.0 mm	330 mm	
1520-0090	1525-0090	9.0 mm	12.7 mm	340 mm	
1520-0095	1525-0095	9.5 mm	13.3 mm	350 mm	
1520-0100	1525-0100	10.0 mm	14.0 mm	360 mm	

Endotracheal Tube, LC Type (Low pressure cuff, Curved):

LC-Magill REF. No.	LC-Murphy REF. No.	Size/I.D.	O.D.	Length	Description
1530-0050	1535-0050	5.0 mm	7.3 mm	240 mm	
1530-0055	1535-0055	5.5 mm	8.0 mm	270 mm	
1530-0060	1535-0060	6.0 mm	8.7 mm	280 mm	*LC-Magill series: Magill tip
1530-0065	1535-0065	6.5 mm	9.3 mm	290 mm	*LC-Murphy series: Murphy eye
1530-0070	1535-0070	7.0 mm	10.0 mm	300 mm	
1530-0075	1535-0075	7.5 mm	10.7 mm	310 mm	Low pressure cuff Curved Graduation Pilot balloon X-ray opaque line
1530-0080	1535-0080	8.0 mm	11.3 mm	320 mm	
1530-0085	1535-0085	8.5 mm	12.0 mm	330 mm	
1530-0090	1535-0090	9.0 mm	12.7 mm	340 mm	
1530-0095	1535-0095	9.5 mm	13.3 mm	350 mm	
1530-0100	1535-0100	10.0 mm	14.0 mm	360 mm	



1520-OXXX

1525-OXXX

1530-OXXX

1535-OXXX

SILICONE ENDOTRACHEAL TUBE







Endotracheal Tube, OW Type (Ordinary cuff, Wire reinforced):

OW-Magill REF. No	OW-Murphy REF. No.	Size/I.D.	O.D.	Length	Description
1540-0045	1545-0045	4.5 mm	6.7 mm	220 mm	
1540-0050	1545-0050	5.0 mm	7.3 mm	240 mm	
1540-0055	1545-0055	5.5 mm	8.0 mm	270 mm	
1540-0060	1545-0060	6.0 mm	8.7 mm	280 mm	*OW-Magill series: Magill tip
1540-0065	1545-0065	6.5 mm	9.3 mm	290 mm	*OW-Murphy series: Murphy eye
1540-0070	1545-0070	7.0 mm	10.0 mm	300 mm	
1540-0075	1545-0075	7.5 mm	10.7 mm	310 mm	Ordinary cuff
1540-0080	1545-0080	8.0 mm	11.3 mm	320 mm	Wire reinforced
1540-0085	1545-0085	8.5 mm	12.0 mm	330 mm	Graduation
1540-0090	1545-0090	9.0 mm	12.7 mm	340 mm	Pilot balloon
1540-0095	1545-0095	9.5 mm	13.3 mm	350 mm	
1540-0100	1545-0100	10.0 mm	14.0 mm	360 mm	

Endotracheal Tube, LW Type (Low pressure cuff, Wire reinforced):

LW-Magill REF. No	LW-Murphy REF. No.	Size/I.D.	O.D.	Length	Description
1550-0050	1555-0050	5.0 mm	7.3 mm	240 mm	
1550-0055	1555-0055	5.5 mm	8.0 mm	270 mm	
1550-0060	1555-0060	6.0 mm	8.7 mm	280 mm	*LW-Magill series: Magill tip
1550-0065	1555-0065	6.5 mm	9.3 mm	290 mm	*LW-Murphy series: Murphy eye
1550-0070	1555-0070	7.0 mm	10.0 mm	300 mm	
1550-0075	1555-0075	7.5 mm	10.7 mm	310 mm	Low pressure cuff
1550-0080	1555-0080	8.0 mm	11.3 mm	320 mm	Wire reinforced
1550-0085	1555-0085	8.5 mm	12.0 mm	330 mm	Graduation
1550-0090	1555-0090	9.0 mm	12.7 mm	340 mm	Pilot balloon
1550-0095	1555-0095	9.5 mm	13.3 mm	350 mm	
1550-0100	1555-0100	10.0 mm	14.0 mm	360 mm	

O: ordinary cuff, L: low pressure cuff, W: wire reinforced



1540-OXXX:LW-MAGILL

1545-OXXX:LW-MYRPHY

1550-OXXX:LW-MAGILL

1555-OXXX:LW-MYRPHY



**fortune
medical**

6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com

