



SILICONE STOMACH (GASTRIC) TUBE

Silicone for Enhanced Patient Comfort




GASTROENTEROLOGIC
SERIES

SILICONE STOMACH (GASTRIC) TUBE

Fortune® Silicone Stomach (Gastric) Tube is designed for nutrition supplement to stomach and may be recommended for various purposes: for patients who cannot take food or swallow, take enough food by mouth to keep nutrition, congenital defects of mouth, esophagus, or stomach.

Fortune® Silicone Stomach (Gastric) Tube Features:

- Medical grade silicone for superior biocompatibility.
- Proper softness and excellent flexibility for enhanced patient comfort.
- Tapered tip for easy insertion.
- Depth marking for position confirmation.




SPECIFICATIONS

Stomach (Gastric) Tube:

REF. No.	Size	Length	Description
2020-0004	4FR		
2020-0005	5FR	400 mm	- Attached plug
2020-0006	6 FR		- Depth marking
2020-0008	8 FR	800 mm	- X-ray opaque line
2020-0010	10 FR		

REF. No.	Size	Length	Description
2020-0012	12 FR		
2020-0014	14 FR		- Attached plug
2020-0016	16 FR	1200 mm	- Depth marking
2020-0018	18 FR		- X-ray opaque line
2020-0020	20 FR		




2020-00XX

NASOGASTRIC TUBE WITH ENFIT CONNECTOR

Fortune® Silicone Nasogastric Tube with ENFit

Connector:


- The silicone nasogastric tube with purple ENFit connector is soft and flexible, It includes a soft tapered tip for easy and safe insertion, and a silicone holder. These components are medical grade silicone for enhancing patient comfort and superior biocompatibility. It can reduce the discomfort contact with the tube.
- Transparent tube for easy visual inspection, X-ray-opaque line for easy and accurate confirmation of tube position by X-ray.
- Depth markings for convenient confirmation of the tube position.
- The nasogastric tube has assembly the ENFit male connector and ENFit female cap when the feeding is stopped. The ENFit female cap can be close by rotated clockwise to prevent the air into the stomach.
- The ENFit connector can prevent the misconnection risk, and provide a safer feeding and medication channel for the patient.



SPECIFICATIONS

Silicone Nasogastric Tube with ENFit Connector:

REF. No.	Size.	Length	Description
2020-0506	6FR	400 mm	
2020-0508	8FR	800 mm	- ENFit Male Connector
2020-0510	10FR		- ENFit Female Cap
2020-0512	12FR		- Size Marking
2020-0514	14FR		- Depth Marking
2020-0516	16FR	1200 mm	- X-Ray Opaque Line
2020-0518	18FR		
2020-0520	20FR		



2020-05XX



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
 Danshuei Dist, New Taipei City. 251 Taiwan
 Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
 www.fortunemed.com Mail:info@fortunemed.com