



SILICONE SUMP DRAIN TUBE



SILICONE SUMP DRAIN TUBE

Fortune® Silicone Sump Drain Tube Features:

- Medical grade silicone for superior compatibility.
- One or two lumens facilitate air circulation and maintain equilibrium of air pressure.
- Each lumen is connected with a filter to protect patient against possible infection from the air when suction and irrigation is proceeding.
- Does not create excessive suction and small capillary to bleed within the wound.
- Proper hardness maintains rigidity for mass suction and irrigation.
- Transparent tube allows easy visual inspection and fluid observation.
- X-ray opaque line.
- Sterilized double packaging.



SPECIFICATIONS


Sump Drain Tube:

REF. No.	Size (H x W)	Length	Eye pattern	Filters
2016-0114	7 x 13 mm	14"	Oval	1
2016-0120	7 x 13 mm	20"	Oval	1
2016-0214	9 x 17 mm	14"	Circular	2
2016-0220	9 x 17 mm	20"	Circular	2
2016-0314	9 x 17 mm	14"	Oval	2
2016-0320	9 x 17 mm	20"	Oval	2


REF. No.	Size (H x W)	Length	Eye pattern	Filters
2016-1314	9 x 17 mm	14"	Oval	0
2016-1320	9 x 17 mm	20"	Oval	0




2016-13XX



2016-01XX



2016-02XX



2016-03XX