



# SILICONE SURGICAL RULER



# SILICONE SURGICAL RULER

## INTENDED USE:

The “Fortune” Silicone Surgical Ruler is used by the surgeon to obtain an approximate measurement of the circumference of organ (e.g., urethra) and the length of organ or tissue in intraoperative.


## FEATURES:

- 100% medical grade silicone for superior biocompatibility.
- Solid silicone strap with markers from 2 cm to 12 cm in 0.5 cm increments.
- Soft silicone material for better comfort.
- Flexible measuring instrument.
- Single use.
- Sterile package.
- Not made with natural rubber latex.



## SPECIFICATIONS

Product	REF. No.	Size.	Description
Silicone Surgical Ruler	1860-7510	Width: 20.5mm Length: 200 mm	Markers from 2 cm to 12 cm in 0.5 cm increments



1860-7510



**fortune**  
medical

6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,  
Danshuei Dist, New Taipei City. 251 Taiwan  
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266  
www.fortunemed.com Mail:info@fortunemed.com

