



SILICONE SURGICAL RULER



SILICONE SURGICAL RULER

INTENDED USE:

The “Fortune®” Silicone Surgical Ruler is used by the surgeon to obtain an approximate measurement of the circumference of organ (e.g., urethra) and the length of organ or tissue in intraoperative.


FEATURES:

- Medical grade silicone for superior biocompatibility.
- Solid silicone strap with markers from 2 cm to 12 cm in 0.5 cm increments.
- Soft silicone material for better comfort.
- Flexible measuring instrument.
- Single use.
- Sterile package.
- Not made with natural rubber latex.



SPECIFICATIONS

Product	REF. No.	Size.	Description
Silicone Surgical Ruler	1860-7510	Width: 20.5mm Length: 200 mm	Markers from 2 cm to 12 cm in 0.5 cm increments



1860-7510



fortune
medical

6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com

