



SILICONE VESSEL ID LOOPS




SILICONE VESSEL ID LOOPS

INTENDED USE:

The “Fortune” Silicone Vessel ID Loops are used in complicated surgical procedures by offering retraction, occlusion, and identification of arteries, nerves, veins, tendons or ureters.

Fortune® Silicone Vessel ID Loops features:

- 100% medical grade silicone for superior biocompatibility.
- Soft and smooth material for better comfort
- Silicone radiopaque
- Non-adherent and non-absorbent to tissue
- Non-toxic and latex-free
- Sterile package
- Available in four sizes and four colors




SPECIFICATIONS

REF. No.	Size.	Application
9005-0503R		Red Color
9005-0503B	Width: 0.5 mm Thickness: 0.3 mm	Blue Color
9005-0503Y	Length: 450 mm	Yellow Color
9005-0503W		White Color

REF. No.	Size.	Application
9005-1309R		Red Color
9005-1309B	Width: 1.3 mm Thickness: 0.9 mm	Blue Color
9005-1309Y	Length: 450 mm	Yellow Color
9005-1309W		White Color



9005-0503X



9005-1309X


GENERAL SURGICAL SERIES

SILICONE VESSEL ID LOOPS

REF. No.	Size.	Application	REF. No.	Size.	Application
9005-2510R		Red Color	9005-5013R		Red Color
9005-2510B	Width: 2.5 mm Thickness: 1.0 mm	Blue Color	9005-5013B	Width: 5.0 mm Thickness: 1.3 mm	Blue Color
9005-2510Y	Length: 450 mm	Yellow Color	9005-5013Y	Length: 450 mm	Yellow Color
9005-2510W		White Color	9005-5013W		White Color



9005-2510X



9005-5013X



fortune
medical

6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com