

INSTRUMENT MAT (NON-STERILE)

The instrument mat can be a neutral zone
which is easy to move.

NON-MAGNETIC INSTRUMENT MAT (NON-STERILE)


Fortune® Non-magnetic Instrument Mat:

Prevent surgical instruments from falling and causing contamination or damage to affect the operation.

Before use, please thoroughly clean, sterilize, and dry the instrument mat.

Feature:


- Medical-grade silicone.
- The pouch on both sides of the non-magnetic Instrument Mat - Pouch Type (9006-0001-02) can hold surgical instruments.
- Surgical instruments are kept within immediate reach of the surgeon's hand and are safely and quickly to retrieve.



SPECIFICATIONS

Non-magnetic Instrument Mat- Sheet Type


REF. No.	Size	Description
9006-0001-01	L540 x W300 mm	- Sheet Type - MR Safe



9006-0001-01

Non-magnetic Instrument Mat- Pouch Type

REF. No.	Size	Description
9006-0001-02	L410 x W300 mm	- Pouch Type - MR Safe



9006-0001-02

MAGNETIC INSTRUMENT MAT (NON-STERILE)

Fortune® Magnetic Instrument Mat:

Prevent surgical instruments from falling and causing contamination or damage to affect the operation. Before use, please thoroughly clean, sterilize, and dry the instrument mat.

Feature:


- Medical-grade silicone.
- It contains magnets that the magnetic attraction provides a more stable retaining force on sloping contoured surfaces.
- Surgical instruments are kept within immediate reach of the surgeon's hand and are safely and quickly to retrieve.
- The surface of the product is smooth, and it is easy to clean.
- It can be rolled and placed in a surgical tray for sterilization with other surgical instruments and reused.



SPECIFICATIONS

Magnetic Instrument Mat

REF. No.	Size	Description
9006-0101	L380 x W300 mm	- MR Unsafe



9006-0101