



SURGICAL INSTRUMENT TRANSFER TRAY (NON-STERILE)



SURGICAL
INSTRUMENT
SERIES


SURGICAL INSTRUMENT TRANSFER TRAY (NON-STERILE)

Fortune® Surgical Instrument Transfer Tray:

The surgical instrument transfer tray will minimize the risk of sharp injuries in operating theatres. Surgical instrument transfer tray acts as a clearly defined neutral zone area for the placement of sharp instruments to transfer and pass sharp instruments by the Hands Free Technique (HFT) during surgery.


Feature:

- Medical-grade silicone.
- Equipped with a handle designed to prevent fingers from entering the tray during handling, ensuring user safety.
- Two sets of instrument holders positioned in opposite directions for optimal convenience in placing and transferring instruments.
- Smooth and Easy-to-Clean Surface.
- Incorporates numerous holes in the tray to promote quick drying, preventing water accumulation.



SPECIFICATIONS

REF. No.	Size	Description
9006-0002	128 x 230 x 24 mm	- Non-sterile



9006-0002



6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
 Danshuei Dist, New Taipei City. 251 Taiwan
 Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
 www.fortunemed.com Mail:info@fortunemed.com