



SILICONE NEPHROSTOMY CATHETER

The Right Choice for Quality Catheters




UROLOGICAL
SERIES

SILICONE NEPHROSTOMY CATHETER

Fortune® Silicone Nephrostomy Catheter Features:

- Medical grade silicone for superior biocompatibility.
- Transparent tube with X-ray opaque line. It allows for easy visual inspection and fluid observation. It also allows for confirmation of intubated tube using X-ray.
- Soft and uniformly inflated balloon makes the tube sit well against the kidney.
- Smooth round shaft can minimize trauma during insertion and withdrawal.
- Four types available for different applications.
- Sterilized double packaging



SPECIFICATIONS

Nephrostomy Catheter, X-ray Opaque Tip, 2 Sided Holes:

Ordinary Balloon REF. No.	Integrated Balloon REF. No.	Size	Balloon	Length	Description
1851-0208	--	8 FR	1.5 cc/ml		
1851-0210	--	10 FR			
1851-0212	4851-0212	12 FR	2 cc/ml		
1851-0214	4851-0214	14 FR	3 cc/ml		- X-ray opaque tip - 2 sided holes
1851-0216	4851-0216	16 FR		420 mm	- Depth marking
1851-0218	4851-0218	18 FR			- X-ray opaque line
1851-0220	4851-0220	20 FR	5 cc/ml		
1851-0222	4851-0222	22 FR			
1851-0224	4851-0224	24 FR			

Nephrostomy Catheter, Open Tip:

Ordinary Balloon REF. No.	Integrated Balloon REF. No.	Size	Balloon	Length	Description
1852-0208	--	8 FR	1.5 cc/ml		
1852-0210	--	10 FR			
1852-0212	4852-0212	12 FR	2 cc/ml		
1852-0214	4852-0214	14 FR	3 cc/ml		- Open tip
1852-0216	4852-0216	16 FR		420 mm	- Depth marking
1852-0218	4852-0218	18 FR			- X-ray opaque line
1852-0220	4852-0220	20 FR	5 cc/ml		
1852-0222	4852-0222	22 FR			
1852-0224	4852-0224	24 FR			



1851-02XX



1852-02XX

SILICONE NEPHROSTOMY CATHETER

Nephrostomy Catheter, Open Tip, 2 Sided Holes:

Ordinary Balloon REF. No.	Integrated Balloon REF. No.	Size	Balloon	Length	Description
1853-0208	--	8 FR	1.5 cc/ml		
1853-0210	--	10 FR			
1853-0212	4853-0212	12 FR	2 cc/ml		
1853-0214	4853-0214	14 FR	3 cc/ml		- Open tip - 2 sided holes
1853-0216	4853-0216	16 FR		330 mm	- Depth marking - X-ray opaque line
1853-0218	4853-0218	18 FR			
1853-0220	4853-0220	20 FR	5 cc/ml		
1853-0222	4853-0222	22 FR			
1853-0224	4853-0224	24 FR			

Nephrostomy Catheter, Open Tip, 2 Sided Holes:

Ordinary Balloon REF. No.	Integrated Balloon REF. No.	Size	Balloon	Length	Description
1854-0208	--	8 FR	1.5 cc/ml		
1854-0210	--	10 FR			
1854-0212	4854-0212	12 FR	2 cc/ml		
1854-0214	4854-0214	14 FR	3 cc/ml		- Open tip - 2 sided holes
1854-0216	4854-0216	16 FR		420 mm	- Depth marking - X-ray opaque line
1854-0218	4854-0218	18 FR			
1854-0220	4854-0220	20 FR	5 cc/ml		
1854-0222	4854-0222	22 FR			
1854-0224	4854-0224	24 FR			



1853-02XX



1854-02XX



fortune
medical

6 F., No. 29, Sec. 2, Jhongjheng E. Rd.,
Danshuei Dist, New Taipei City. 251 Taiwan
Tel.:(886) 2-2624-2233 Fax.:(886) 2-2624-2266
www.fortunemed.com Mail:info@fortunemed.com

